

# MULTITOX ASSAY

## REDUCTION OF DRUG ATTRITION AND ZEBRAFISH

The pharmaceutical industry is under growing pressure to improve R&D productivity to sustain sufficient innovation to replace the loss of revenues due to patent expiration for successful products. Combined with the low number of new molecular entities (NMEs) entering into clinical phases, **new screening strategies** are demanded. Compounds fail for many reasons, but some are more avoidable such as poor oral bioavailability, pharmacokinetic properties or toxicity, and low margins of safety.

A new strategy to reduce attrition and maintain the number of NMEs entering clinical phases requires improving toxicity characterization at early stages of drug discovery combined with more efficient and cost-effective assays.

The **zebrafish** is emerging as a complement to existing *in vitro* technologies and established preclinical *in vivo* models that can be scaled for high-throughput. Technological innovation has helped the zebrafish embryo gain ground as a disease model and an assay system for drug screening. Moreover, the zebrafish embryo offers a cost effective system that combines many features that are desirable for the development of new approaches to drug development.

By using Zebrafish toxicity assays not only for predicting toxicities in later stages but applying them at the very early stages of the Drug Discovery process, it would be feasible to improve the selection of safe candidates or to understand potential toxicities that should be evaluated, thereby decreasing the attrition of drug candidates in more advanced phases due to toxicity.

## MULTITOXICITY ASSAY

Based on the organ attrition related to specific toxicities and its incidence in the Drug Development process, we propose an innovative approach sequentially combining different toxicity assays, to rapidly address the deselection of compounds through toxicity assessment in zebrafish embryos. This approach will also offer the reduction of costs linked to a reduction of compounds tested at each sequential assay.

Biobide has developed a multi-assay that combines three mayor toxicity assays responsible for attrition in the R&D process: Cardiotoxicity, Neurotoxicity and Teratogenesis, followed by a refinement step that is open to the specific requirements of each screening program (Figure 1).




Figure 1. Scheme of the assays and conditions used for the MultiToxAssay

## VALIDATION RESULTS

Out of 56 tested reference compounds, 9 compounds did not show any toxicity at the end of the process (Table 1). In Figure 2 the number of compounds selected in each step of the assay are shown. The assay has an overall sensitivity of 90% and specificity of 100%.

Table 1. List of reference compounds used in the validation process

Classification of products	Reference compounds	Classification of productos	Reference compounds
Toxins	Aflatoxin		Propranolol
Nootropic	Deprenyl	Beta-blockers	Timolol
Analgesic	Acetaminophe		Sotalol
	Lovastatin		Thioridazine
Statins	Simvastatin		Haloperidol
	Mevastatin		Pimozide
		Antipsychotics	Sertindole
Anti-inflammatory/immunosuppressant	Dexamethasone		Risperidone
Antibiotics	Tetracyclin		Carbamazepine
	Penicillin G		Valproic acid
Anticoagulant	Warfarin	Anti-convulsant	Pilocarpine
Gastrointestinal agents	Cisapride		Nicotine
	Verapamil	Colinergic agonist	Tacrine
Calcium blockers	Diltiazem		Astemizole
	Nitrendipine		Terfenadine
	Terodiline	Antihistaminics	Halofantrine
Antitumorals	Tamoxifen		Foscarnet
	Chlorambucil		Mefloquine
	5-Fluorouracil	Anti-infective	Isoniazid
	Thalidomide		Testosterone
	Hydroxyurea		Estradiol
Antidepressants	Fluoxetine	Hormones	Hydrocortisone
	Amitriptyline		Ascorbic acid
Antiarrhythmics	Lidocaine		NAC
	Flecainide	Vitamins/anti-oxidants	Sucrose
	Propafenone		Indirubin-3-oxine
	Amiodarone		Digitoxin
	Disopiramide	Others	
Ina openers	sdz-201106		Ketanserine

Specificity: 90%  
Sensitivity: 100%




Figure 2. Percentage and the number of compounds selected in the multiassay approach for a better selection of candidates.

- ✓ The MultiTox assay can help to decrease the attrition rate linked to toxicity using *in vivo* data at early stages of the Drug Discovery Process by revealing which compounds do not induce any of the main toxicities.
- ✓ The use of this multiassay approach would result in a 50% reduction of the time and cost needed to perform each assay individually.
- ✓ The MultiTox assay is modular and can be customized to any requirement. It allows changing the type, the sequence and even the conditions of each toxicity assay selected for the multiassay approach.